

FRONT DETECTOR ( )

SIGNAL 1  
Data rate: 20 Hz  
Type: test plot  
Save Data: Off  
Zero: 0.0 (Off)  
Range: 0  
Fast Peaks: Off  
Attenuation: 0

BACK DETECTOR ( )

SIGNAL 2  
Data rate: 20 Hz  
Type: test plot  
Save Data: Off  
Zero: 0.0 (Off)  
Range: 0  
Fast Peaks: Off  
Attenuation: 0

COLUMN COMP 1  
(No Detectors Installed)

COLUMN COMP 2  
(No Detectors Installed)

THERMAL AUX 2  
Use: MSD Transfer Line Heater  
Description:  
Initial temp: 300 'C (On)  
Initial time: 0.00 min  
# Rate Final temp Final time  
1 0.0(Off)

POST RUN  
Post Time: 0.00 min

TIME TABLE  
Time

Specifier

Parameter & Setpoint

7673 Injector

Front Injector:  
Sample Washes 0  
Sample Pumps 3  
Injection Volume 1.0 microliters  
Syringe Size 10.0 microliters  
PostInj Solvent A Washes 3  
PostInj Solvent B Washes 3  
Viscosity Delay 3 seconds  
Plunger Speed Fast  
PreInjection Dwell 0.00 minutes  
PostInjection Dwell 0.00 minutes

Back Injector:  
No parameters specified  
Column 1 Inventory Number : ?  
Column 2 Inventory Number :

MS ACQUISITION PARAMETERS

General Information

Tune File : atune.u  
Acquisition Mode : Scan

MS Information

Solvent Delay : 9.00 min  
EM Absolute : False  
EM Offset : 0  
Resulting EM Voltage : 1811.8

[Scan Parameters]

Low Mass : 50.0

Method: MAN\_52.M

Mon Jul 24 14:19:20 2006