

TOPLEVEL PARAMETERS

Method Information For: D:\MSDCHEM\1\METHODS\MAN27.M

Method Sections To Run:

- (X) Save Copy of Method With Data
- () Pre-Run Cmd/Macro =
- (X) Data Acquisition
- (X) Data Analysis
- () Post-Run Cmd/Macro =

Method Comments:

Quantification du rapport Testosterone /Epitestosterone MSD20 injection en split

END OF TOPLEVEL PARAMETERS

INSTRUMENT CONTROL PARAMETERS

Sample Inlet: GC
Injection Source: GC ALS
Mass Spectrometer: Enabled

6890 GC METHOD

OVEN

Initial temp: 160 'C (On) Maximum temp: 325 'C
Initial time: 0.00 min Equilibration time: 0.50 min
Ramps:
Rate Final temp Final time
1 4.00 255 0.00
2 30.00 300 2.75
3 0.0(Off)
Post temp: 0 'C
Post time: 0.00 min
Run time: 28.00 min

FRONT INLET (UNKNOWN)

BACK INLET ()

Mode: Split
Initial temp: 280 'C (On)
Pressure: 175.0 kPa (On)
Split ratio: 12.1:1
Split flow: 11.7 mL/min
Total flow: 15.2 mL/min
Gas saver: Off
Gas type: Helium

COLUMN 1

COLUMN 2

Method: MAN27.M

Sat Jul 22 13:31:13 2006

Page:1

USADA 0201

188