

TOPLEVEL PARAMETERS

Method Information For: C:\MSDCHEM\1\METHODS\MAN_52.M
 Method Sections To Run:

- Save Copy of Method With Data
- MSTOP Pre-Run Cmd/Macro =
- Instrument Control Pre-Run Cmd/Macro =
- Data Analysis Pre-Run Cmd/Macro =
- Data Acquisition
- Data Analysis
- MSTOP Post-Run Cmd/Macro =
- Instrument Control Post-Run Cmd/Macro =
- Data Analysis Post-Run Cmd/Macro =

Method Comments:
 This is the default method

END OF TOPLEVEL PARAMETERS

INSTRUMENT CONTROL PARAMETERS

6890 GC METHOD

OVEN

Initial temp: 70 'C (On)
 Initial time: 1.00 min

Maximum temp: 325 'C
 Equilibration time: 0.50 min

Ramps:

#	Rate	Final temp	Final time
1	30.00	270	12.00
2	10.00	300	3.00
3	0.0(Off)		

Post temp: 70 'C
 Post time: 0.00 min
 Run time: 25.67 min

FRONT INLET (SPLIT/SPLITLESS)

Mode: Splitless
 Initial temp: 280 'C (On)
 Pressure: 145.0 kPa (On)
 Purge flow: 20.0 mL/min
 Purge time: 2.00 min
 Total flow: 25.4 mL/min
 Gas saver: Off
 Gas type: Helium

BACK INLET (UNKNOWN)

COLUMN 1

Capillary Column
 Model Number: Agilent 19091s-433
 Max temperature: 325 'C
 Nominal length: 30.0 m
 Nominal diameter: 250.00 um
 Nominal film thickness: 0.25 um
 Mode: constant pressure
 Pressure: 145.0 kPa
 Nominal initial flow: 2.3 mL/min
 Average velocity: 56 cm/sec
 Inlet: Front Inlet
 Outlet: MSD
 Outlet pressure: vacuum

COLUMN 2
 (not installed)

Method: MAN_52.M

Mon Jul 24 14:19:20 2006

Page: 1