

DATA ANALYSIS PARAMETERS

Method Name: D:\MSDCHEM\1\METHODS\MAN27.M

Percent Report Settings

Sort By: Signal

Output Destination

Screen: No
Printer: Yes
File: No

Integration Events: AutoIntegrate

Generate Report During Run Method: No

Signal Correlation Window: 0.020

Qualitative Report Settings

Peak Location of Unknown: Apex

| Library to Search | Minimum Quality |
|-------------------|-----------------|
| DEMO.L | 0 |

Integration Events: AutoIntegrate

Report Type: Summary

Output Destination

Screen: No
Printer: Yes
File: No

Generate Report During Run Method: No

Quantitative Report Settings

Report Type: Summary

Output Destination

Screen: Yes
Printer: No
File: No

Generate Report During Run Method: No

Method: MAN27.M

Mon Jul 24 12:55:02 2006

Page:1

USADA 0206

191