

TOPLEVEL PARAMETERS  
-----

Method Information For: C:\MSDCHEM\1\METHODS\MAN\_52.M  
Method Sections To Run:

- Save Copy of Method With Data
- MSTOP Pre-Run Cmd/Macro =
- Instrument Control Pre-Run Cmd/Macro =
- Data Analysis Pre-Run Cmd/Macro =
- Data Acquisition
- Data Analysis
- MSTOP Post-Run Cmd/Macro =
- Instrument Control Post-Run Cmd/Macro =
- Data Analysis Post-Run Cmd/Macro =

Method Comments:  
This is the default method

END OF TOPLEVEL PARAMETERS  
-----

INSTRUMENT CONTROL PARAMETERS  
-----

=====

6890 GC METHOD

=====

OVEN

Initial temp: 70 'C (On)  
Initial time: 1.00 min

Maximum temp: 325 'C  
Equilibration time: 0.50 min

Ramps:

#	Rate	Final temp	Final time
1	30.00	270	12.00
2	10.00	300	3.00
3	0.0 (Off)		

Post temp: 70 'C  
Post time: 0.00 min  
Run time: 25.67 min

FRONT INLET (SPLIT/SPLITLESS)

Mode: Splitless  
Initial temp: 280 'C (On)  
Pressure: 145.0 kPa (On)  
Purge flow: 20.0 mL/min  
Purge time: 2.00 min  
Total flow: 25.4 mL/min  
Gas saver: Off  
Gas type: Helium

BACK INLET (UNKNOWN)

COLUMN 1

Capillary Column  
Model Number: Agilent 19091s-433  
Max temperature: 325 'C  
Nominal length: 30.0 m  
Nominal diameter: 250.00 um  
Nominal film thickness: 0.25 um  
Mode: constant pressure  
Pressure: 145.0 kPa  
Nominal initial flow: 2.3 mL/min  
Average velocity: 56 cm/sec  
Inlet: Front Inlet  
Outlet: MSD  
Outlet pressure: vacuum

COLUMN 2  
(not installed)

Method: MAN\_52.M

Fri Aug 04 10:24:41 2006

Page: 1