

High Mass : 550.0
Threshold : 150
Sample # : 2 A/D Samples 4
Plot 2 low mass : 50.0
Plot 2 high mass : 550.0

[MSZones]

MS Quad : 150 C maximum 200 C
MS Source : 230 C maximum 250 C

END OF MS ACQUISITION PARAMETERS

END OF INSTRUMENT CONTROL PARAMETERS

DATA ANALYSIS PARAMETERS

Method Name: C:\MSDCHEM\1\METHODS\MAN_52.M

Percent Report Settings

Sort By: Signal

Output Destination
Screen: No
Printer: Yes
File: No

Integration Events: AutoIntegrate

Generate Report During Run Method: No

Signal Correlation Window: 0.020

Qualitative Report Settings

Peak Location of Unknown: Apex

Library to Search Minimum Quality
C:\temp\IRMS.L 90

Integration Events: AutoIntegrate

Report Type: Summary

Output Destination
Screen: No
Printer: Yes
File: No

Generate Report During Run Method: No

Quantitative Report Settings

Method: MAN_52.M

Fri Aug 04 10:24:41 2006

Page: 3